

Not just customers. Fans!

Featuring patented NASA spacesuit technology, Game Ready is chosen by medical practitioners, professional teams, and athletes for the treatment of a variety of injuries and post-surgery rehabilitation. The Game Ready System is being used by more than 100 pro sports teams, 600 professional athletes, and 200 colleges and universities nationwide.



"Game Ready is my cold/compression device of choice, out of all the products on the market. I've used it consistently, with fabulous results. I recommend it as a device for every trainer to have in their training room."

—Jason Powell
Head ATC
Los Angeles Clippers



"The Game Ready device is excellent...it is the best modality that has come out in quite some time."

—Joe Carroll, PT MS SCS
Sports Clinical Specialist
Owner Cape Cod Rehabilitation




"Game Ready has become an important addition to our patients' post-operative recovery."

—Sanford Kunkel, MD
Orthopedic Surgeon



"The patients are singing its praises. I have spoken to two patients the morning after their surgery and they couldn't say enough."

—Maggie Cladwell, RN
Beltway Surgery Center



**Injury?
Surgery?**
**Get Game Ready.
And get back in
the game.**

Interested in learning more?
For a Game Ready Representative in your area
or for more information, Call toll-free
1.888.Game.Ready (1.888.426.3732)
www.gameready.com

Come back like the Pros!

Are you recovering from an injury or surgery? Do what professional athletes are doing. Get Game Ready. When you're sidelined following an injury or orthopedic surgery, you've got one goal. You don't want the couch. You want to get back in action. You want to get Game Ready.



Get Game Ready

RICE (Rest, Ice, Compression, Elevation) has long been used to treat acute and chronic injury and assist in orthopedic surgery rehabilitation. Game Ready offers adjustable cold therapy and intermittent compression in **one easy-to-use system**, making it easy for you to apply the two most difficult-to-manage aspects of the RICE regimen. Many professional athletes have discovered the Game Ready System and choose it over traditional methods of applying RICE because Game Ready's form-fitting, flexible Wraps make it easy to deliver cold and intermittent compression to most major joints. Get back into the game—with Game Ready.



"Our medical staff has been looking for a high-quality and durable product that provides cold and compression for years, and we finally found it! The results have been outstanding and it has become our modality of choice for acute and chronic injuries."

—Chuck Barta
ATC, Head Athletic Trainer
Minnesota Vikings




The Game Ready Difference.

How Does Game Ready Work?

Game Ready's **anatomic Wraps** are each designed to fit a specific body part: the ankle, knee, back, hip/groin, shoulder, elbow and wrist. The Wrap is applied to the injured area and secured. The **Control Unit** is filled with ice and water and connected to the Wrap. When the unit is turned on, the temperature is adjustable, while microprocessors enable both time and the intermittent compression that is applied to the injured area via the Wrap.

What Makes Game Ready Different?

Want to get over an injury or surgery and back in the game? Whether you are an athlete or active person, surgery patient, or medical professional—you should consider the Game Ready System. Game Ready offers the simultaneous delivery of both **adjustable cold therapy and intermittent compression** in one easy-to-use unit. The portable Game Ready System, with its flexible, circumferential, dual-action Wraps and its adjustable settings of time, temperature and compression, really stands out from the crowd.



Rest. Game Ready. Elevate.

Finally, there's an easy-to-use application of the two most difficult-to-manage aspects of RICE (Rest, Ice, Compression, Elevation). Now it's easy to treat injuries and assist in post-operative rehabilitation by applying the RICE regimen: **Rest, Game Ready, Elevate!**

1} ADJUSTABLE THERAPY

Game Ready's adjustable technology enables continuous, rapid circulation of ice water - along with intermittent compression - throughout circumferential Wraps.

2} ANATOMIC WRAPS

Game Ready does not take a one-wrap-fits-all approach. Our anatomically designed Dual-Action Wraps are built to treat specific parts of the body for a close, comfortable fit.

3} EASE OF USE

Elastic bandages and messy ice packs just can't compete with Game Ready's adjustable system and no-hassle wraps. Easy to apply and comfortable to use, athletes and patients often ask for Game Ready by name.

4} REIMBURSEMENT

May be available through health insurance carriers for the use of Game Ready. For more information or assistance with reimbursement, please contact us toll-free at 1.888.426.3732.

6} PORTABILITY

Game Ready's Carry Bag and Wrap Bag help you to conveniently and safely store and transport your important therapy equipment.

7} ONE YEAR WARRANTY

One year warranty on the Control Unit and the Wrap's air and water chamber.



Built for comfort and convenience.

Where elastic bandages, messy ice packs, and other alternatives offer limited adjustability, the Game Ready System allows you to easily adjust pressure, manage temperature, and select treatment time settings—with just the push of a few buttons and the twist of a dial. Now you can simply and simultaneously apply intermittent compression and adjustable cold therapy to deliver the two most difficult-to-manage aspects of the RICE regimen in one easy-to-use treatment system.



Five Easy Steps



1

Fill the Control Unit with ice and water, then connect to a power source.



2

Apply the Wrap to the injured area and secure.



3

Attach the Connector Hose first to the Wrap, then to the Control Unit.



4a

If you have a Pre-Programmed Control Unit, press Power, adjust to 3 Snowflakes, then press the Timer button to select the desired program.



4b


If you have a standard Control Unit, press Power, adjust to 3 Snowflakes, then select pressure and time settings.




5

Press Start.

Spacesuit technology. Dual-action wraps.



Game Ready Two-Part Wrap



Game Ready Wrap Technology

Flexible Wrap Technology

Developed by scientists and doctors using **NASA spacesuit technology**, Game Ready's flexible wraps conform to the body for circumferential cold and compression therapy. Designed for most major body parts, the portable Game Ready System featuring **patented form-fitting, no-hassle wraps**, can assist in your injury treatment or post-surgery rehabilitation.

Dual-Action Process

The Game Ready Wrap is designed in two parts, featuring a removable outer sleeve and an inner **dual-action chamber**. The dual-action air and water chamber enables water to flow rapidly through the Wrap, while the air inflates and deflates to deliver active cold and compression therapy.

 "This system is worth its weight in gold!"
—Kenny Roberts, Jr
Professional Motorcycle Champion

Circumferential application, anatomic design.

Game Ready's flexible fabric Wraps surround the affected area for good coverage and a customized fit. Unlike many products that use a one-design-fits-all approach, the Game Ready System features a comprehensive range of light-weight, easy-to-use Wraps that are specifically engineered for all major body parts and available in a variety of sizes.

ankle wrap
(2 sizes)



knee wrap

shoulder wrap
R & L, (2 sizes)



elbow wrap

hip/groin wrap
(R & L)



back wrap

wrist wrap



articulated knee wrap